

## Appendix 11

### Recommendations for size of bronchoscope with different airway devices

Airway device size	New scope sizes
ETT 3.5	<2.7
ETT 5.5	2.7
ETT 6.5	2.7
ETT 8.0	5.3
ETT 9.0	7.3
ETT 10.5	8.7

ETT – Endotracheal Tube

**Table adapted from:**

*Brimacombe J, Dunbar-Reid K. The effect of introducing fiberoptic bronchoscopes on gas flow in laryngeal masks and tracheal tubes. Anaesthesia. 1996 Oct;51(10):923-8*